

# (PRACTICAL) COMPUTATIONAL PHYSICS

Physics 551  
Lecture 14

Extra Reading

Optional Exercise

Recommended

# NOTATION

- This lecture slides for this course will attempt to use a uniform notation throughout. A normal paragraph looks like this.
- ☞ *Italicized paragraphs with pen bullets will indicate definitions, with the defined word or phrase shown in **SMALL-CAPS**.*
- ☞ Pencil bullets will indicate the introduction of **new notation**.
- ☞ Pointing hand bullets indicate important points that might otherwise be overlooked.

# ANNOUNCEMENTS

# ANNOUNCEMENTS

- To **clone** this week's **LabVIEW demonstration materials** please open the ***Git Shell*** utility and invoke

```
$ git clone https://github.com/hughdickinson/CompPhysL14Labview.git
C:/Users/ComputationalPhysics/Documents/labview/lecture14
```

- ☞ You can also find this command on the Blackboard Learn website.

# ANNOUNCEMENTS

- **Project Deadline:**
  - **Midnight May 4th 2015** (although you may want to finish earlier if you have exams).
- **Project Submission:**
  - Create a **new directory** for your project called *FinalProject* in your private homework repositories on **GitHub**.
  - Add **all** of the elements of your project (including the manual) to this directory.
  - I will clone the contents of this directory **immediately** after the submission deadline.

# ANNOUNCEMENTS

- **Project Requirements:**
  - **Functional** executable programs for C++ or LabVIEW.
  - **Functional** Python Scripts and Modules.
  - **Fully documented** and **annotated** code (Doxygen-compatible for C++).
  - A **user manual** including instructions for program **invocation** and **usage** and a **simple tutorial** section.


# ANNOUNCEMENTS

- **Project Tips:**
  - Document code **as you write it.**
  - It is always easier to write documentation when the while the program logic is clear in your mind!
  - Use **git** to **regularly commit** your code as you develop.
  - Protect your work by regularly pushing committed changes to GitHub.
  - You may wish to back up your work in other ways.

# SPECIALIZED LABVIEW DATA TYPES

# PATH DATA

- LabVIEW provides a **generic** *Open/Create/Replace File* VI to **open** files that contain data in an **arbitrary** format.



**Predefined File  
Constants**

# PATH DATA

- The *OCR File* VI accepts an optional input connection that has the special **path** data type.
  - ☞ If an explicit path is not connected to the appropriate input then the *OCR File* VI will present a **file browser dialog** to the user.
- The *OCR File* VI returns a **File Reference** that can be connected to other LabVIEW file I/O VIs.
  - ☞ LabVIEW defines several *constant path values* that can be used to **assemble** file paths programmatically.

# DEMONSTRATION



## Opening Files and Assembling Paths Using Predefined Constants

# WAVEFORM DATA

- Many LabVIEW applications acquire data from instruments at a **constant rate**.
  - ☞ In LabVIEW such data are typically stored as **WAVEFORM** data types.
  - ☞ Waveforms are aggregate data types that comprise **at least three** distinct components.
    - An **array of sample values**.
    - A value corresponding to the **time of the first sample**.
    - A value specifying the **time between consecutive samples**.

# WAVEFORM DATA

- Tabview provides VIs to associate arbitrary **attributes** with waveform data.
- **Specific** attributes can be used to influence how LabVIEW **plotting indicators** render **waveform** data.




# DEMONSTRATION

Working with **Waveform Data**

# DYNAMIC DATA

- ☞ The LabVIEW package provides several, highlight configurable **EXPRESS VIs** that typically output **DYNAMIC TYPE DATA**.
- ☞ LabVIEW provides VIs to **explicitly convert** dynamic data to fundamental data types.
- ☞ Many of the **plotting indicators** provided by LabVIEW accept dynamic data as an input - no explicit conversion is required.
- Dynamic data types can introduce **computational overhead**, and are typically **only** be used in conjunction with **Express VIs**.
- Like **Waveform** data, dynamic data types can also have **arbitrary attributes** associated with them.


# DYNAMIC DATA



Again, **specific** attributes can be used to influence how LabVIEW **plotting indicators** render **dynamic** data.

**Explicit  
Conversion VIs**

**Attribute  
Manipulation VIs**




# DEMONSTRATION

Working with **Dynamic Data**

# VARIANT DATA

- ☞ The **VARIANT** data type is designed to introduce limited type-flexibility with LabVIEW's strict typing model.
- Variant data stores a **value** and **metadata** specifying how the value should be interpreted.
- ☞ VIs exist to **safely convert** variant data to arbitrary fundamental types. Impossible conversions produce **errors**, not improper values.
- “Flattening” VIs exist that permit the data and metadata that comprise a variant to be examined.

# VARIANT DATA




# DEMONSTRATION

Working with **Variant Data**

# ERROR CLUSTERS

- Error handling in LabVIEW involves **special** instances of the **cluster** data types.

- ☞ LabVIEW **ERROR CLUSTERS** comprise three **specific** elements
  - A **boolean** value with identifier **status**.
  - A **32 bit integer** value with identifier **code**.
  - A **string** value with identifier **source**.




# HANDLING ERRORS IN LABVIEW


# LABVIEW ERROR HANDLING

- If any of the sub-VIs or functions that populate the block diagram generate an error, the **default** behaviour of LabVIEW entails **three** operations.
  - Suspension of program execution.
  - Highlighting the point in the program where the error occurred.
  - Displaying an error dialog box.
- ☞ This default behaviour can be **disabled** via the LabVIEW user interface.


# LABVIEW ERROR HANDLING



# LABVIEW ERROR HANDLING



# LABVIEW ERROR HANDLING



- After **disabling** the default behaviour, the **Dialog & User Interface** palette provides VIs to implement **customized** error-handling mechanisms.


**Error Dialog VIs**

**Error Cluster  
Creation VIs**

**Error Manipulation  
VIs**

# LABVIEW ERROR HANDLING

- A **subset** of LabVIEW VIs report errors using bare **numeric codes** instead of Error Clusters.
- The *Error Cluster From Error Code VI* builds **Error Clusters** from **Error Codes**.
- The error code is used to populate the **code** element of the **Error Cluster**.
- Multiple Error Clusters can be **compounded** using the extensible *Merge Errors* function.



# LABVIEW ERROR HANDLING

- Error Clusters can be **defined** using context-specific **integer** Error Codes.
- LabVIEW **reserves** specific integer ranges for internal use.
- There are **three allowed ranges** for developer-defined Error Codes

-8999

→

-8000

8000

→


8999

500000

→

599999

## Error Definition VIs



# DEMONSTRATION

## Error Handling in LabVIEW

# FILE INPUT AND OUTPUT

# FILE INPUT AND OUTPUT

- ☞ LabVIEW provides **several** VIs that are tailored to **read** and **write** data in numerous **proprietary** and **generic** formats.
  - All of these VIs are accessible from the **File I/O** palette of the **block diagram** context menu.
  - VIs exist to handle **arbitrary textual** and **binary** data as well as formatted textual data including **XML** and **spreadsheets**.
- ☞ **Dedicated** VIs are provided for reading and writing **Waveform** data.
- ☞ Support is also provided for **file compression** using the **Zip** protocol.

# FILE INPUT AND OUTPUT

**Spreadsheet**  
File Handling

**Generic Text**  
File Handling

**Proprietary**  
**LabVIEW**  
**TDMS**


**Waveform-Specific**  
File Handling

**Formatted Text**  
File Handling

**Generic Binary**  
File Handling

**XML**  
File Handling

**Zip**  
File Handling



# DEMONSTRATION

## **File Input** and **Output** in LabVIEW

# PLOTTING DATA

# PLOTTING USING WAVEFORM CHARTS

- ☞ LabVIEW provides **CHART** indicators to display data **as it is accumulated**.
- **Charts** plot and maintain values in a **fixed-length** buffer that stores a **limited** number of the most recently provided data.
- The data are managed in a **first-in, first-out** (FIFO) **queue**. When the data buffer is full, the **oldest data are discarded first**.
- **Charts** are useful for **monitoring** continuous data acquisition.

# PLOTTING USING WAVEFORM CHARTS


- Several **numeric** types and various **layouts** of data can be connected to a **Chart** input.
  - ☞ The **layout** of the data defines how the **values** should be **interpreted** and **plotted**.
- The **simplest** possible way to provide data to a **Chart** is to input a **single** numeric datum **whenever an update is required**.
- Alternatively, multiple values can be supplied in a single update by connecting an **array** of numeric types to the **Chart** input.

# PLOTTING USING WAVEFORM CHARTS

- **Waveform** data can also be connected to LabVIEW **Chart**-type indicators.
  - ☞ The **start timestamp** and **sampling interval** provided by the waveform data type are used to properly calibrate the **Chart**'s x-axis.
  - ☞ If **successive** waveforms are not strictly **temporally adjacent** then only the most recent waveform is plotted.
- Otherwise, successive wavelengths are **appended** to the plot.

# PLOTTING USING WAVEFORM CHARTS

- Single numeric data, arrays and waveform data are rendered as single traces by LabVIEW charts.
- LabVIEW clusters can be used to provide data to charts that is rendered as multiple traces.



# DEMONSTRATION

## **Plotting** using **Charts** in LabVIEW

# PLOTTING USING WAVEFORM GRAPHS

- ☞ LabVIEW provides **GRAPH** indicators to display data **after** it is accumulated.
- **Graphs** do **not** maintain a history of previous updates. Only the **most recent** update is displayed.
- **Graphs** are useful for **plotting** summaries of accumulated data.
- **Graphs** allow the **x-coordinates** of plotted data to be explicitly specified.
- LabVIEW provides several different types of **Graph** indicator.


# PLOTTING USING WAVEFORM GRAPHS

- ☞ LabVIEW provides **WAVEFORM GRAPH** indicators to display **evenly sampled** data e.g. a snapshot from a **Waveform Chart**.
- **Waveform Graph** indicators only plot single-valued functions of the x-coordinate i.e.  $y = f(x)$ .
- ☞ **Waveform Graph** indicators do not require **explicit specification of x-coordinate values**. They will be inferred from the **types** and **layout** of the data that are supplied.
- ☞ **Waveform Graphs** also use data layout to determine whether the data values should be rendered as **single** or **multiple** traces.



# PLOTTING USING WAVEFORM GRAPHS

- The illustrated data layouts (involving combinations of **single data, arrays, clusters** and **waveforms**) can be used to plot **single** traces, **without** explicitly specifying the x-coordinates.

X-values set to  
**array index**




X-values computed using  
**start values** and **sampling intervals**




# PLOTTING USING WAVEFORM GRAPHS

- The illustrated data layouts (involving combinations of **single data**, **arrays**, **clusters** and **waveforms**) can be used to plot **multiple** traces that share **identical** x-coordinates.

X-values set to  
**innermost**  
array index




X-values computed using start value ( $x_0$ )  
and sampling interval ( $\delta x$ ).  
All traces share the **same** ( $x_0, \delta x$ ).



# PLOTTING USING WAVEFORM GRAPHS


- The illustrated data layouts (involving combinations of **single data, arrays, clusters** and **waveforms**) can be used to plot **multiple** traces that have **different** x-coordinates or different numbers of points.

X-values set to  
**innermost**  
array index




X-values computed using start value ( **$x_0$** )  
and sampling interval ( **$\delta x$** ).

*Same Sampling*



*Different Sampling*



# DEMONSTRATION

## **Plotting using Waveform Graphs in LabVIEW**

# PLOTTING USING INTENSITY CHARTS AND GRAPHS

- ☞ LabVIEW provides **INTENSITY CHART** indicators to display **three-dimensional** data on a **two-dimensional** Cartesian plane **as they are accumulated.**
- ☞ LabVIEW provides **INTENSITY GRAPH** indicators to display **three-dimensional** data on a **two-dimensional** Cartesian plane **after they are accumulated.**
- 👉 **Intensity Charts** and **Intensity Graphs** accept a **two-dimensional array** of numeric data as input.

# PLOTTING USING INTENSITY CHARTS

- The **indices** of the input data array are interpreted as **x** and **y** coordinates, implying **even sampling** in those dimensions.
  - ☞ By **default**, the **column** indices of the input array are interpreted as **y**-coordinates and the **row** indices are interpreted as **x**-coordinates.
- The **values** of the array elements are mapped to a color value and used to render a **colored region** at the corresponding x and y coordinates.

# PLOTTING USING INTENSITY CHARTS

- ☞ **Intensity Charts** behave like **Waveform Charts** and maintain a **limited** history of recently provided data.
- Each **new** three-dimensional data set is rendered at the **rightmost** edge of the plot and earlier data are shifted to the left.
- **Intensity Graphs** handle input data in a similar manner to **Waveform Graphs**.
- ☞ If new data are supplied to an **Intensity Graph**, previously supplied data are **discarded**.

# DEMONSTRATION

## **Plotting using Intensity Charts and Intensity Graphs in LabVIEW**

# PLOTTING USING XY GRAPHS




- ☞ LabVIEW provides **XY GRAPHS** to render **arbitrarily sampled two-dimensional** data sets in a **Cartesian plane**.
- 👉 XY Graphs **require explicit** specification of the **x** and **y** coordinates for each datum that they render.
- 👉 Several **numeric** types and various **layouts** of data can be connected to an XY Graphs as input.
- 👉 The **layout** of the data defines how the **values** should be **interpreted** and **plotted**.

# PLOTTING USING XY GRAPHS

- ☞ LabVIEW provides **XY GRAPHS** to render **arbitrarily sampled two-dimensional** data sets in a **Cartesian plane**.
- 👉 XY Graphs **require explicit** specification of the **x** and **y** coordinates for each datum that they render.
- 👉 Several **numeric** types and various **layouts** of data can be connected to an XY Graphs as input.
- 👉 The **layout** of the data defines how the **values** should be **interpreted** and whether they should be **plotted** as single or multiple data sets.

# PLOTTING USING XY GRAPHS


- The illustrated data layouts (involving combinations of **single numeric data**, **complex valued data**, **arrays** and **clusters**) can be supplied to **XY Graphs** to plot **single** datasets.
- When supplying data as **complex values**, the **real** part is interpreted as the **x**-coordinate and the **imaginary** part is interpreted as the **y**-coordinate.





# PLOTTING USING XY GRAPHS

- The illustrated data layouts (involving combinations of **single numeric data**, **complex valued data**, **arrays** and **clusters**) can be supplied to **XY Graphs** to plot **multiple** datasets.

*Dataset 1*



*Dataset 2*



# DEMONSTRATION

## **Plotting** using **XY Graphs** in LabVIEW


# PLOTTING USING 2D AND 3D GRAPHS

- ☞ To handle alternative rendering options for multidimensional data, LabVIEW defines **abstract 2D GRAPH** and **3D GRAPH** concepts.
  - **Concrete implementations** of the 2D Graph and 3D Graph concepts render 2-dimensional and 3-dimensional data in **specific** ways.
- ☞ Each concrete implementation provides a **HELPER** VIs that **format input data** appropriately.
- ☞ LabVIEW implements **graphs with error bars** as concrete implementations of the 2D Graph concept.

# PLOTTING USING 2D AND 3D GRAPHS

*Data  
Generation*

*Concrete 2D Graph  
Implementation Plots Data*



*Helper VI  
Formats Data*

# DEMONSTRATION

## **Plotting** using **2D** and **3D Graphs** in LabVIEW

# RECOMMENDED READING

## **Learn LabVIEW Online Video Tutorials**

(<http://www.ni.com/academic/students/learn-labview/>)

## **LabVIEW Basics Online Reference**

(<http://www.ni.com/getting-started/labview-basics/>)

## **LabVIEW Manual: Types of Graphs and Charts**

([http://zone.ni.com/reference/en-XX/help/371361L-01/lvconcepts/types\\_of\\_graphs\\_and\\_charts/](http://zone.ni.com/reference/en-XX/help/371361L-01/lvconcepts/types_of_graphs_and_charts/))

## **LabVIEW Manual: Handling Errors**

([http://zone.ni.com/reference/en-XX/help/371361L-01/lvconcepts/error\\_checking\\_and\\_error\\_handling/](http://zone.ni.com/reference/en-XX/help/371361L-01/lvconcepts/error_checking_and_error_handling/))

## **LabVIEW Manual: Using the Dynamic Data Type**

([http://zone.ni.com/reference/en-XX/help/371361L-01/lvconcepts/dynamic\\_data\\_type/](http://zone.ni.com/reference/en-XX/help/371361L-01/lvconcepts/dynamic_data_type/))

# LECTURE 14 HOMEWORK

Review the **Recommended Reading** items listed on the previous slide.

Continue to refine your **final project proposal** and begin working on your **final project** once it is approved.